

Application Brief

cMET/EGFR Copy Number Variation and cMET Gene Expression Panel

Unique All-In-One-Well Copy Number Variation and Gene Expression Multimodal Assay


INTRODUCTION

cMET is a proto-oncogene that encodes for Hepatocyte Growth Factor Receptor (HGFR), a receptor tyrosine kinase, which plays an essential role in normal cellular function and oncogenesis. Recent studies have identified cMET as a biomarker for various cancers as well as for drug resistance in case of anti-EGFR therapy. cMET can be shown to highly over-express and show significant copy number variations under such situations. Additionally, copy number variation in EGFR can occur by a complex interaction between cMET and EGFR. With so many biomarkers in so many complex pathways targeted by so many drugs, there is a need for all-in-one-well, multiplex and multimodal tests. Using the ICEPlex® system, we created an 18-target, multimodal quantitative cMET/EGFR Copy Number Variation and cMET Gene Expression Panel designed to detect:

- cMET/EGFR gene copy number variation (CNV)
- Chromosomal copy number variation
- Over-expression of cMET mRNA

SUMMARY

- Delivers accurate and precise quantification of copy number variation and gene expression.
- Provides built-in controls such as calibration and process controls to assure quality of the materials and samples.
- Enables amplification from FFPE tissue materials because of the smaller amplicon (<100 nucleotides) and requires minimum nucleic acid input.
- Tests versatile sample types: FFPE blocks, fresh-frozen tissues, or cell lines.



UNIQUE MULTIMODAL CAPABILITY

TYPICAL DATA

A real-time, all-in-one-well 18-target assay was developed to determine cMET/EGFR CNV and cMET expression on the ICEPlex system. Performance testing using total nucleic acids from cell lines with known cMET CNV and overexpression as well as chromosome 7 polysomy matched published results. Resulting data were similar for fresh-frozen or FFPE-derived material. Low amounts of FFPE-extracted material were sufficient to accurately determine cMET/EGFR CNV and cMET expression status.

NOTE: Data were normalized to reference genes Ref2 and Ref3 and fold-change calculated.

Figure 1




Figure 2




Figure 1 Computer-generated gel view showing simultaneous detection of cMET/EGFR and cMET expression targets.

Figure 2. cMET/EGFR copy number variation, expression and polysomy, all-in-one-well.

METHOD HIGHLIGHTS

- This 18-target panel was constructed to quantify cMET/EGFR Gene Copy Number and cMET Gene Expression.
- For cMET expression profiling, three primer pairs amplifying three different regions of cMET mRNA were designed.
- For cMET/EGFR CNV three primer pairs each amplifying three different genomic regions of cMET and EGFR were designed.
- For normalization and relative quantification, two reference genes were used; two primer pairs for each reference gene mRNA were designed for cMET mRNA quantification and one primer pair for each reference gene was designed for cMET/EGFR copy number variation.
- Multiplexed reactions were optimized on the ICEPlex system using total nucleic acids from various cell lines.
- The fluorescently labeled amplicons for the different cMET/EGFR CNV and cMET gene expression profile targets were injected, separated and detected in the capillary electrophoresis module of the ICEPlex system.
- Amplification curves for all targets and controls were generated automatically by the ICEPlex software and cycle thresholds (Cts) were calculated.

FOR MORE INFORMATION

For a list of publications and to find out more about how PrimeraDx can help your lab, please contact us at 508.618.2300 or visit www.primeradx.com.

The ICEPlex system and the ICEPlex cMET/EGFR Copy Number Variation and cMET Gene Expression Assay are for Research Use Only and have not been approved for in vitro diagnostic use by the FDA. The presented information is for demonstration purposes only.